

# Table A-26

Quality Control Acceptance Criteria for Method E418.1 — TPH

|   |  |
|---|--|
| <b>Water</b><br><b>RL-Water</b> = 0.05 mg/L<br><b>Accuracy Water (% R)</b> = 80-120<br><b>Precision Water (RPD)</b> = ±20 | <b>Soil</b><br><b>RL-Soil</b> = 3.2 mg/kg<br><b>Accuracy Soil (%R)</b> = 75-125<br><b>Precision Soil (RPD)</b> = ±20 |
|---|--|

| QC Check                                      | Minimum Frequency  | Acceptance Criteria  | Corrective Action <sup>a</sup>  | Flagging Criteria <sup>b</sup>   |
|---|--|--|---|--|
| Initial calibration (5 standards and a blank) | Daily (prior to sample analysis)   | Correlation coefficient (r) > 0.995<br>Calibration MUST meet acceptance criteria prior to sample analysis. | Identify and repeat outlying point(s); recalculate curve using repeated point(s).                                 | Apply R to all results for specific analyte(s) for all samples associated with the calibration   |
| ICV   | Daily following initial calibration  | ± 10% true value   | 1) Repeat ICV.<br><br>2) If not compliant, evaluate system and correct problem.<br>3) Repeat initial calibration. | Apply J to positive results and UJ to non-detects all results for specific analyte(s) for all samples associated with the calibration      |
| Continuing calibration verification (CCV)     | Daily, before sample analysis, after every 10 samples, and at the end of each batch. | + 10% true value   | 1) Identify and correct problem.<br><br>2) Recalibrate and reanalyze all samples since last valid CCV.            | Apply J to positive results and UJ for non-detects for all results for specific analyte(s) for all samples associated with the calibration |

| QC Check                        | Minimum Frequency                                | Acceptance Criteria | Corrective Action <sup>a</sup>  | Flagging Criteria <sup>b</sup>  |
|---------------------------------|--|---------------------|---|---|
| Method blank                    | 1 per analytical batch and per preparation batch | < RL                | 1) If sample analyte concentration is < RL or if the sample analyte concentration is > 10 times the concentration in the method blank, then report results and annotate in case narrative   | Apply U to all results for the specific analyte in all samples in the associated analytical batch if concentration is less than 5 times blank concentration.  |
|                                 |  |                     | 2) If preparative method blank does not meet item 1), re extract/re-analyze if still within HT and enough sample volume; if not within HT or enough sample, contact project chemist for decision.   |   |
| Laboratory control sample (LCS) | 1 per preparation batch and analytical batch     | 80-120%             | 1) If the preparative LCS recovers high outside the acceptance criteria and the analyte is ND, flag the LCS results and annotate in the case narrative.<br>2) If the preparative LCS fails the acceptance criteria (other than shown in item 1), re-extraction and reanalysis will be necessary if samples are still within holding time and enough sample volume; if not, contact the project chemist for a decision | For specific analyte in all samples in the associated analytical batch: if the LCS %R > UCL, apply J to all positive results<br><br>if the LCS %R < LCL, apply J to all positive results, apply UJ to all non-detects. If LCS <50%, R flag results. |

| QC Check                                     | Minimum Frequency                | Acceptance Criteria  | Corrective Action <sup>a</sup>   | Flagging Criteria <sup>b</sup>  |
|--|----------------------------------|--|--|---|
| Matrix spike/matrix spike duplicate (MS/MSD) | One MS/MSD per every 40 samples. | 75-125% recovery and RPD < 20%                                     | If the MS and/or MSD is outside of either accuracy or precision tolerances and LCS results are acceptable, flag MS/MSD results and annotate in the case narrative. | For specific analyte in all samples in the associated analytical batch: if the MS/MSD is > UCL, J all hits; if < LCL but greater than 30%, J all hits, UJ all non-detects. If < 30%, R all non-detects and J all hits. If precision outside criteria, J all hits, UJ all non-detects. |
| MDL study                                    | Once per 12 month period         | Detection limits established shall be < 1/2 the RLs in this table. | None   | Apply R to all results for specific analyte(s) in all samples analyzed.   |
| Results reported between MDL and RL          | None                             | None   | None   | Apply J to all results between MDL and RL   |

<sup>a</sup> All corrective actions associated with project work shall be documented, and all records shall be maintained by the laboratory.

<sup>b</sup> Flagging criteria are applied when acceptance criteria were not met and corrective action was not successful or corrective action was not performed.